

Sulfo-Cy5 dUTP

Chemical Properties

CAS No. :

Formula: C₅₀H₆₁Li₄N₆O₂₂P₃S₂

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cy5 dUTP is a modified triphosphate nucleotide utilized for DNase labeling. This compound serves as a bright, far-red emitting dye specifically designed for the Cy5 channel.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481