

Sulfo-Cy3.5 amine

Chemical Properties

CAS No. :

Formula: C44H52K2N4O13S4

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cy3.5 amine is a sulfonated version of Cy3.5, characterized by excellent water solubility. Its excitation and emission maxima are at 591 nm and 604 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481