

Cy3.5 alkyne

Chemical Properties

CAS No. :

Formula: C₄₁H₄₄ClN₃O

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3.5 alkyne is a fluorescent dye and a click chemistry reagent containing an alkyne group. It can undergo a copper-catalyzed azide-alkyne cycloaddition (CuAAC) with azide-containing molecules. This compound is useful for labeling and detecting azide-containing molecules/biomolecules via fluorescence spectroscopy following the cycloaddition reaction.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481