

Sulfo-Cy7.5 DBCO

Chemical Properties

CAS No. :

Formula: C₆₁H₅₇K₃N₄O₁₄S₄

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sulfo-Cy7.5 DBCO is a dye derivative of Cyanine 7.5 (Cy7.5) that contains a DBCO group. The presence of sulfonate ions enhances its water solubility, making it suitable for aqueous applications. Sulfo-Cy7.5 DBCO can bind to biomolecules such as proteins and antibodies to track their location and dynamics in biological samples.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481