

## Cy5.5 phosphoramidite

### Chemical Properties

CAS No. :

Formula: C<sub>55</sub>H<sub>71</sub>ClN<sub>5</sub>O<sub>3</sub>P

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cy5.5 phosphoramidite, a cyanine derivative, serves as a fluorescent labeling reagent employed in the synthesis of fluorescently labeled oligonucleotides.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481