

ICG-TCO

## Chemical Properties

CAS No. :

Formula: C<sub>56</sub>H<sub>68</sub>N<sub>4</sub>O<sub>6</sub>S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ICG-TCO is a fluorescent dye composed of indocyanine green (ICG) conjugated with TCO. It can be utilized for the detection and visualization of molecules containing tetrazines.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481