

JNK-IN-15, Cell-Permeable

Chemical Properties

CAS No. :

Formula: C191H336N70O48S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	JNK-IN-15, Cell-Permeable (JNK inhibitor III), is an inhibitor of JNK. It is applicable in the research of age-related neurodegenerative diseases.
Targets(IC50)	JNK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481