

JNK-IN-15, Cell-Permeable, Negative Control

Chemical Properties

CAS No. :

Formula: C191H336N70O48S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | JNK-IN-15, Cell-Permeable, Negative Control serves as the negative control for JNK-IN-15, Cell-Permeable. JNK-IN-15, Cell-Permeable functions as an inhibitor of JNK. |
| Targets(IC50) | JNK |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481