

$\alpha$ -Secretase Substrate II, Fluorogenic

## Chemical Properties

CAS No. :

Formula: C100H144N30O20S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	$\alpha$ -Secretase Substrate II, Fluorogenic is an internally quenched fluorescent peptide substrate for $\alpha$ -secretase, featuring the $\alpha$ -secretase cleavage site of $\beta$ -amyloid precursor protein (APP). Its excitation/emission (Ex/Em) wavelengths are 340/490 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481