

11R-CaN-CON

Chemical Properties

CAS No. :

Formula: C135H248N70O34S2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	The compound 11R-CaN-CON is an inactive control peptide used as a negative control for 11R-CaN-AID.
Targets(IC50)	Phosphatase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481