

Pen-N3

Chemical Properties

CAS No. :

Formula: C51H78N10O15

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Pen-N3 is a PDZ peptide that inhibits Wnt/b-catenin signaling by recognizing disordered proteins.
Targets(IC50)	Wnt/beta-catenin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481