

Arg-Arg-Arg-Arg acetate

Chemical Properties

CAS No. :

Formula: C₂₄H₅₀N₁₆O₅.xC₂H₄O₂

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Arg-Arg-Arg-Arg (acetate), composed of four arginine units, is applicable in the development of gene delivery vectors that utilize cell-penetrating peptides.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481