

## E3 Ligase Ligand-linker Conjugate 92

## Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>33</sub>N<sub>5</sub>O<sub>6</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	E3LigaseLigand-linker Conjugate 92 is a conjugate of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), composed of Thalidomide and the corresponding linker. E3LigaseLigand-linker Conjugate 92 functions as a Cereblon ligand to recruit the CRBN protein and serves as a crucial intermediate in the synthesis of complete PROTAC molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481