

## E3 Ligase Ligand-linker Conjugate 99

### Chemical Properties

CAS No. :

Formula: C<sub>28</sub>H<sub>37</sub>N<sub>5</sub>O<sub>6</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	E3LigaseLigand-linker Conjugate 99 is a conjugate of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), composed of Thalidomide and its corresponding linker. It functions as a Cereblon ligand to recruit CRBN protein, and serves as a crucial intermediate in synthesizing complete PROTACs molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481