

E3 Ligase Ligand-linker Conjugate 100

Chemical Properties

CAS No. :

Formula: C₂₆H₃₁N₅O₆

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3LigaseLigand-linker Conjugate 100 is a conjugate of an E3 ligase ligand with a linker (E3LigaseLigand-Linker Conjugates), composed of Thalidomide and its corresponding linker. E3LigaseLigand-linker Conjugate 100 functions as a Cereblon ligand to recruit CRBN protein and serves as a crucial intermediate in the synthesis of complete PROTAC molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481