

## E3 Ligase Ligand-linker Conjugate 103

## Chemical Properties

CAS No. :

Formula: C40H59N7O6S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	E3LigaseLigand-linker Conjugate 103 is a conjugate of an E3 ligase ligand and a linker, comprising (S,R,S)-AHPC and its corresponding linker. (S,R,S)-AHPC is a von Hippel-Lindau (VHL) ligand based on VH032, capable of recruiting the VHL protein. E3LigaseLigand-linker Conjugate 103 serves as a crucial intermediate for the synthesis of complete PROTAC molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481