

## E3 Ligase Ligand-linker Conjugate 104

### Chemical Properties

CAS No. :

Formula: C33H45N5O7

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	E3LigaseLigand-linker Conjugate 104 is a combination of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), consisting of Thalidomide and the associated linker. This conjugate functions as a Cereblon ligand to recruit CRBN protein and serves as a crucial intermediate for the synthesis of complete PROTACs molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481