

PROTAC BTK Degradar-8

Chemical Properties

CAS No. :

Formula: C80H94F2N14O20P2

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PROTACBTK Degradar-8 (compound 3) is an efficient PROTAC BTK degrader with a linker that can couple with ADC antibodies to form PAC.
Targets(IC50)	BTK,PROTAC-Linker Conjugates for PAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481