

TRX-COBI

## Chemical Properties

CAS No. :

Formula: C38H43F3IN3O7

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TRX-COBI ((R,R)-TRX-COBI; TRX-cobimetinib) is an MEK inhibitor that targets the TRX fragment based on Fe <sup>2+</sup> . It exhibits antitumor activity.
Targets(IC50)	MEK

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481