

PBX-7011 mesylate

## Chemical Properties

CAS No. :

Formula: C<sub>25</sub>H<sub>25</sub>N<sub>3</sub>O<sub>9</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PBX-7011 mesylate is an active derivative of camptothecin (Camptothecin). It binds to the DDX5 protein within cells, leading to cell death through the degradation of the DDX5 protein.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481