

SrtA-IN-2

Chemical Properties

CAS No. :

Formula: C₂₀H₁₅FN₂O₄

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Antibacterialagent 175 (compound Y40) is an ML346 analogue and an inhibitor of sortase A (SrtA). It demonstrates inhibitory activity against Staphylococcus aureus SrtA and shows inhibition of biofilm formation. Antibacterialagent 175 acts as an antitoxin against S. aureus infections.
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481