

## EGFR/STAT3-IN-1

### Chemical Properties

CAS No. :

Formula: C<sub>30</sub>H<sub>24</sub>ClFN<sub>4</sub>O<sub>2</sub>S<sub>2</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	EGFR/STAT3-IN-1 (Compound 5k) is a dual inhibitor targeting both EGFR and STAT3, with IC <sub>50</sub> values of 41 nM and 3.34 nM, respectively. It exhibits antitumor activity.
Targets(IC <sub>50</sub> )	EGFR,STAT

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481