


Proroxan hydrochloride

Chemical Properties

CAS No. :	33025-33-1
Formula:	C ₂₁ H ₂₄ ClNO ₃
Molecular Weight:	373.87
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Proroxan hydrochloride is a non-selective adrenergic α -antagonist. Proroxan hydrochloride is utilized in hypertension research.
Targets(IC50)	Adrenergic Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6747 mL	13.3736 mL	26.7473 mL
5 mM	0.5349 mL	2.6747 mL	5.3495 mL
10 mM	0.2675 mL	1.3374 mL	2.6747 mL
50 mM	0.0535 mL	0.2675 mL	0.5349 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481