

Lactoperoxidase, bovine milk

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lactoperoxidase is an antibacterial protein from the mammalian heme peroxidase family. It catalyzes the oxidation of its preferred halide substrates, such as I ⁻ , Br ⁻ , F ⁻ , or the pseudohalide SCN ⁻ , using H ₂ O ₂ .
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481