

TAMRA-isoADPr

Chemical Properties

CAS No. :

Formula: C₅₁H₆₂N₁₂O₁₈P₂

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAMRA-isoADPr is a fluorescently labeled iso-ADP-ribose tracer compound. It serves as a high-affinity probe for the RNF146 WWE domain, with a dissociation constant (K _d) of 45.2 nM.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481