

JJC12-009

## Chemical Properties

CAS No. :

Formula: C<sub>29</sub>H<sub>35</sub>ClN<sub>2</sub>O<sub>2</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	JJC12-009 is a dopamine transporter (DAT) inhibitor, with a K <sub>i</sub> of 5.45 nM.
Targets(IC <sub>50</sub> )	Dopamine Receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481