

DCLK1-IN-3

## Chemical Properties

CAS No. :

Formula: C<sub>23</sub>H<sub>22</sub>N<sub>6</sub>O<sub>5</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DCLK1-IN-3 is a potent and selective inhibitor of doublecortin-like kinase 1 (DCLK1) with an IC <sub>50</sub> of 47 nM, and it is applicable in cancer research.
Targets(IC <sub>50</sub> )	Serine Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481