

DCLK1-IN-4

Chemical Properties

CAS No. :

Formula: C₂₄H₂₄N₆O₅

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DCLK1-IN-4 (compound D2) is a potent and selective inhibitor of doublecortin-like kinase 1 (DCLK1), with an IC ₅₀ of 70 nM. It is applicable to cancer research.
Targets(IC ₅₀)	Serine Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481