

## (S)-Thalidomide-piperazine-(1S,4r)-cyclohexane-NH-Boc

## Chemical Properties

CAS No. :

Formula: C<sub>29</sub>H<sub>39</sub>N<sub>5</sub>O<sub>6</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(S)-Thalidomide-piperazine-(1S,4r)-cyclohexane-NH-Boc is a conjugate of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), composed of Thalidomide and the corresponding Linker. It functions as a Cereblon ligand to recruit CRBN protein and serves as a crucial intermediate for synthesizing complete PROTACs molecules.
Targets(IC50)	Autophagy,Ligands for E3 Ligase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481