

PROTAC PTK6 ligand-O-C2-O-piperidine-Boc

Chemical Properties

CAS No. :

Formula: C₃₆H₄₉N₅O₈S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PROTACPTK6 ligand-O-C2-O-piperidine-Boc is a E3 Ligase Ligand-Linker Conjugates). It consists of (S,R,S)-AHPC and the corresponding linker. The (S,R,S)-AHPC is a von Hippel-Lindau (VHL) ligand based on VH032, which can recruit VHL proteins. This compound serves as a crucial intermediate for the synthesis of complete PROTACs molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481