

E3 ligase Ligand-Linker Conjugate 31

Chemical Properties

CAS No. :

Formula: C43H65N7O7S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3ligaseLigand-Linker Conjugate 31 is a conjugate of an E3 ubiquitin ligase ligand and a linker, which can be utilized in the synthesis of complete PROTACs molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481