

E3 ligase Ligand-Linker Conjugate 32

Chemical Properties

CAS No. :

Formula: C41H63N7O6S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3ligaseLigand-Linker Conjugate 32 is a conjugate of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), composed of (S,R,S)-AHPC and its corresponding Linker. (S,R,S)-AHPC is a von Hippel-Lindau (VHL) ligand based on VH032, capable of recruiting VHL protein. E3ligaseLigand-Linker Conjugate 32 serves as a crucial intermediate in the synthesis of complete PROTAC molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481