

## E3 ligase Ligand-Linker Conjugate 38

## Chemical Properties

CAS No. :

Formula: C<sub>22</sub>H<sub>27</sub>N<sub>5</sub>O<sub>5</sub>

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	E3ligaseLigand-Linker Conjugate 38 is a conjugate of an E3 ligase ligand and a linker, comprising Thalidomide and the corresponding linker. E3ligaseLigand-Linker Conjugate 38 serves as a Cereblon ligand to recruit CRBN proteins and acts as a crucial intermediate in the synthesis of complete PROTACs molecules.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481