

NFL-NH2

## Chemical Properties

CAS No. :

Formula: C31H32ClN3O

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NFL-NH2 is a near-infrared ratiometric fluorescent probe targeting mitochondria. It enables rapid detection of NO levels related to the degree of inflammatory damage in rheumatoid arthritis (RA) mouse models through ratiometric fluorescence imaging. The excitation and emission wavelengths are 650 nm and 780 nm, respectively.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481