

PDI-IN-2

Chemical Properties

CAS No. :

Formula: C₂₄H₂₄O₁₃

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|----------------------------|--|
| Description | PDI-IN-2 (Compound 22) is an inhibitor of PDI or recombinant PDI, with IC ₅₀ values of 0.62 and 0.63 μM, respectively. It also exhibits antitumor activity. |
| Targets(IC ₅₀) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481