

pFBC

Chemical Properties

CAS No. :

Formula: C₁₁H₁₀FN₃O

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	pFBC ([¹⁸ F]pEBC) is a covalent CLIP-tag radioactive tracer utilized for identifying viral reporter gene transfer in the brains of mice, facilitating research in neurobiology.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481