

SIRT4-IN-1

Chemical Properties

CAS No. :

Formula: C19H13N5O6S3

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SIRT4-IN-1 (compound 69) is a selective inhibitor of Sirtuin 4 (Sirt4) with an IC50 value of 16 µM. It does not exhibit significant effects on other sirtuin isoforms.
Targets(IC50)	Sirtuin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481