

## 3-Tosylimidazolidine-2,4-dione

### Chemical Properties

CAS No. :

Formula: C<sub>10</sub>H<sub>10</sub>N<sub>2</sub>O<sub>4</sub>S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	3-Tosylimidazolidine-2,4-dione (compound 2C) is a hydantoin derivative that exhibits inhibitory activity against the Pf3D7 strain, with an IC <sub>50</sub> value of 3.97 nM.
Targets(IC <sub>50</sub> )	Parasite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481