

FAP Ligand 2

Chemical Properties

CAS No. :

Formula: C43H46F2N11NaO7S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | FAP Ligand 2 is a derivative of FAPI-46 and serves as a specific tracer targeting fibroblast activation protein (FAP). It exhibits excellent tumor uptake, making it a promising molecular tracer for FAP tumor imaging. |
| Targets(IC50) | Others |
| In vitro | FAP Ligand 2 exhibits a high affinity for FAP. |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481