

Cy5.5(Me)-C3-DBCO

Chemical Properties

CAS No. :

Formula: C₅₈H₅₇ClN₄O₂

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy5.5(Me)-C3-DBCO is a click chemistry reagent containing a cycloalkyne group. The DBCO moiety facilitates copper-free bioorthogonal click chemistry, characterized by rapid reaction kinetics and excellent stability.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481