

BDP TMR ceramide

Chemical Properties

CAS No. :

Formula: C₃₉H₅₆BF₂N₃O₄

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BDP TMR ceramide is a ceramide labeled with BDP TMR, a water-soluble fluorescent dye with high lipophilicity (Ex: 542 nm, Em: 572 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481