

phospho-STAT3-IN-2

Chemical Properties

CAS No. :

Formula: C18H15NO3

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Phospho-STAT3-IN-2 (compound 4D) is a STAT3 inhibitor that effectively inhibits the phosphorylation of STAT3. In mouse xenograft tumor models, phospho-STAT3-IN-2 significantly reduces tumor size without exhibiting drug toxicity to other organs or tissues.
Targets(IC50)	STAT

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481