

N3-PC

## Chemical Properties

CAS No. :

Formula: C43H83N4O8P

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	N3-PC (Azido-PC) is a photoactivatable phosphatidylcholine used in the synthesis of cell membrane structures. It is a click chemistry reagent featuring an azido group.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481