

TAAR1 agonist 1

Chemical Properties

CAS No. :

Formula: C₁₃H₁₄N₂O₂

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TAAR1 agonist 1 (compound 6E) is a potent agonist of TAAR1, demonstrating effective activation of both the TAAR1-Gs and TAAR1-Gq pathways. It significantly alleviates schizophrenia-like cognitive deficits induced by MK-801.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481