

KRAS ligand 4

Chemical Properties

CAS No. :

Formula: C₂₄H₂₈ClF₃N₆O₃

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	KRAS ligand 4 (Compound 2) is a dual-functional molecular glue degrader based on SOS1. It reduces downstream signaling markers pERK and pS6, and exhibits antiproliferative activity in cells with various KRAS mutations.
Targets(IC50)	Molecular Glues,PROTACs,Ras

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481