

## Gαq/11 protein-IN-1

### Chemical Properties

CAS No. :

Formula: C<sub>19</sub>H<sub>27</sub>N<sub>5</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Gαq/11 protein-IN-1 (compound F33) is a Gαq/11 protein inhibitor with an IC <sub>50</sub> value of 9.4 μM, exhibiting anti-tumor activity.
Targets(IC <sub>50</sub> )	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481