

MNK-IN-4

Chemical Properties

CAS No. :

Formula: C₁₅H₁₇N₅

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MNK-IN-4 (compound D25) is a potent and selective MNK inhibitor used in treating acute spleen injury related to sepsis.
Targets(IC50)	MNK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481