

PROTAC SOS1 degrader-6

Chemical Properties

CAS No. : 3032282-51-9


Formula: C46H55F3N8O6S

Molecular Weight: 905.04

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PROTACSOS1 degrader-6 (compound 23) is an effective degrader of SOS1 PROTACs and works synergistically with KRASG12C inhibitors.
Targets(IC50)	PROTACs,Ras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1049 mL	5.5246 mL	11.0492 mL
5 mM	0.221 mL	1.1049 mL	2.2098 mL
10 mM	0.1105 mL	0.5525 mL	1.1049 mL
50 mM	0.0221 mL	0.1105 mL	0.221 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481