

PROTAC PAPD5 degrader 1

Chemical Properties

CAS No. :

Formula: C49H63N5O16

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PROTACPAPD5 degrader 1 (compound 12b) demonstrates the ability to inhibit hepatitis A virus (HAV) and hepatitis B virus (HBV) both in vitro and in vivo, with an IC50 of 10.59 μ M and a CC50 greater than 50 μ M in Huh7 cells.
Targets(IC50)	HBV,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481