

PD-1/PD-L1-IN-42

Chemical Properties

CAS No. :

Formula: C85H129N21O24S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PD-1/PD-L1-IN-42 (Compound B8.4) is an inhibitor of the PD-1/PD-L1 interaction, with an EC50 of 0.1 μ M. It is applicable in cancer immunotherapy research.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481